

Installing MagAO-X On The Telescope

This procedure describes how to install the MagAO-X instrument on the Magellan Clay Telescope

Initial Conditions


- [] Instrument in LCO cleanroom on legs

Preparations

- [] Shutdown MagAO-X
- [] Turn off table legs and remove air connection
- [] Remove all cables
- See detailed procedure for removing 2K DM cables
- [] Electronics Rack
- [] Ensure that roll-out shelves are restrained
- [] If not installed, Install side panels
- [] Close and lock doors
- [] Tape keys down
- [] Instrument
- [] Remove eyepiece
- [] Turn off blower, and remove hose
- [] Tape over any exposed holes (from cables, etc)
- [] Secure any loose cables
- [] Shrink wrap the instrument
- [] Install solar blanked over shrink wrap
- [] Cart and Rigging
- [] Verify all cart hardware is in-hand
- [] Verify two wire harnesses are in hand
- [] Partially assemble cart, leaving one long side off
- [] Store card out of the way so table can roll into unpacking area

Rig Onto Cart

- [] Lower table legs onto the casters by turning the 16 leveling bolts, and remove the metal pads
- [] Roll the instrument out of the cleanroom
- [] At least 4 people required, 2 to push, 2 to hold the cleanroom sides open
- [] Continue rolling the instrument through the garage door into the unpacking area
- [] Move the cart, currently with 3 sides, around the legs
- NOTE: be careful to not bump the legs with the cart
- [] Attach the 4th side of the cart
- [] Ensure that the 8 large bolts on the cart are snug but not tight
- [] Attach the lifting wire-harness to each side of the cart
- [] Attach the load spreader with straight extensions to the crane, using a crane scale



The load spreader attached to the crane for lifting the cart

- [] Place the load spreader in the center position (the cart is symmetric)
- [] Lift the load spreader, and position it over the instrument

- [] Being careful to not bump the instrument, lower the load spreader and attach the lifting harness D rings. Use 2x shackles to extend the length to reach the cart on the floor.



Lifting harnesses attached with shackle extensions

- [] Position a person at each corner of the cart

NOTE: Do not allow the cart to bump the legs or the table uncontrolled

- [] Slowly lift the cart (**320 lbs**) until it is touching the bottom of the table



The cart being lifted to the bottom of the table.

- [] Install the 8 bolts attaching the cart to the table. Loosen bolts on the cart as needed.
- [] Once the cart is bolted to the table bottom, while **320 lbs** is still on the crane, tighten all cart bolts. Do not over-tighten, make 1/4 turn after the washers are no longer free. This is to avoid excessive stress on the table.
- [] Install the triangle stabilizing ropes between the crane hook and the lifting fixture IAW the below figure. Tighten, but do not cause them to pick the load.



*The triangle stabilizing ropes should be tight, but not become the lifting point for the load.

CAUTION: be sure that the load spread does not hit the cart when being repositioned. - ☐ Reposition the load spreader center to the instrument + cart position marked on it.

- ☐ Ensure that there is room to move the legs out from under the table, opening the garage door into the cleanroom if necessary.
- ☐ Position a person at each corner of the cart to stabilize it during the lift.
- ☐ Position two people to remove the legs from under the table
- ☐ Lift the table off the legs. The weight of the table + cart is **1920 lbs**
- ☐ Move the legs out from under the table.



The cart and instrument ready to be set down on the wheels, with legs out of the way.

- [] Set the cart down on its wheels.

Transport MagAO-X To Clay



MagAO-X will be strapped to the Isuzu similarly to how the ASM is transported.



MagAO-X will be unloaded at the telescope using the lift gate, adjusted for slow operation as it is for the ASM.

Install MagAO-X On The Platform

Transport Electronics

Install Electronics and Cable

Install AOC in Control Room